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EN INSTRUCTIONS FOR USE

DESCRIPTION

The ATS Bulb Set is a complete single-use postoperative surgical drainage system for collection and recovery of autologous blood for immediate autotransfusion. The system is available in two versions with different levels of negative suction pressure.

- The ATS Bulb Set has a 50 mmHg (6,6 KPa) initial pressure and an average operating pressure of about 30 mmHg (4,0 KPa);
- The ATS Bulb Set with inox inner spring has a 90 mmHg (12,0 KPa) initial pressure and an average operating pressure of about 65 mmHg (8,6 KPa).

Moreover, an Autotransfusion Unit (bag) is available separately; for its use see the section "Autotransfusion Unit replacement" in these directions for use.

The system is composed of the following parts (see Figure 1):

- A silicone bulb with a 400 ml capacity that is highly biocompatible and transparent to have a clear vision of drained blood;
- A 200 µm filter situated inside the silicone bulb to filter the macroaggregates;
- Double no return valve: one situated at the bulb entrance and the other at the Autotransfusion Unit entrance;
- A unidirectional valve, which is normally closed, to drain the air from the bulb in order to keep the latter pressed, maintaining the vacuum activated;
- Autotransfusion Unit with a 40 µm filter inside that allows eliminating microaggregates before the reinfusion;
- Optional components for the completion of the system: one or more drainage tubes, with or without trocar, an extension connecting the bulb to the bag.

Note: This device, and its components, can only be used in safety conditions in the field of applications and using the procedures indicated on this operating instruction sheet for the specific type of product. The manufacturer declines all liability relating to improper use or use different to that indicated. This instruction leaflet must accompany the device during its entire lifetime and must always be at hand for consultation.

GENERAL WARNINGS

Do not use if the package has been opened or damaged. Carefully read this instruction leaflet before using the product. Single-use product. Reuse may lead to alteration of performance and risks of cross-contamination.

CONTRAINDICATIONS

Blood recovery and reinfusion are contraindicated in the following cases: anomalous hepatic and/or renal functionality; malignant wounds; contamination and/or sepsis; use of fluids unsuitable for reinfusion (ex. Betadine, etc.); Presence of amniotic fluid or bile; presence of haemostatic agents, coagulation problems. It is unadvisable the reinfusion of the recovered blood beyond 6 hours from the onset of the recovery (ref.Guidelines AABB). The infusion of the recovered blood is also contraindicated in all those medical situations considered inadequate by the physician in charge. The above mentioned contraindications refer to the immediate autotransfusion of autologous blood and are independent of the Ats Bulb Set device.

PREPARING THE SYSTEM

The ATS Bulb Set system and its components are packed in a double sterile wrapping. Open the package using aseptic procedures.

1. Position the drainage tubes in the wound following an aseptic procedure, proceeding from inside to outside. Fix the drainage tubes by means of suture or plaster.
2. Connect the connecting tube (1) to the patient drainage through the Y connector provided and all possible available connectors. Cut the multi-size connectors in correspondence with the drainage measurement (the measurement is indicated on the connector section); insert the drainage on the connector (see Figure 2).
Note: A multi-size connector is supplied already open, with a CH 6 measure, and it is identified by the relevant "OPEN" label.
3. To put the suction into operation, close the blood recovery clamp (13) and press the silicone bulb (see Figure 3). The air comes out from the unidirectional valve (4) while the no return valve (5) prevents that the former moves to the patient. Insert the cap on the unidirectional valve (4).
4. Apply the Autotransfusion Unit on the way out of the silicone bulb (7). Indicate patient identity on the Autotransfusion Unit and the time of the start of blood collection. As recovered blood is kept at room temperature, the time between the starting of the collection and the starting of the reinfusion must be within and not later than 6 hours.
5. Open the output clamp (13). Once the maximum of 400 ml recovery capacity has been obtained, the direct communication between the silicone bulb and the Autotransfusion Unit allows passing automatically to a mechanic suction phase and to operate by gravity. The drainage continuity is a safety guarantee that reduces the possibility of huge haematoma and the risk for infections. To obtain this operational continuity, the Ats Bulb Set must be positioned below the patient level.

RECOVERY OF AUTOLOGOUS COLLECTED BLOOD

The silicone bulb transparency (6) visualizes the quantity and quality of blood passed in the system through the macroaggregate filter. To drain the collected blood in the bulb, press slowly with both hands until the whole content has been transferred to the Autotransfusion Unit (see Figure 3). Once this operation is over, the bulb maintains a suction condition ready to continue the drainage.

AUTOTRANSFUSION UNIT REPLACEMENT (BAG)



1. Close the clamp (13).
2. Close the clamp (14) of the Autotransfusion Unit.
3. Remove, detaching the Autotransfusion Unit, and close with the cap (15).
4. The ATS Bulb Set uses an Autotransfusion Unit of 1,000 ml capacity that is normally sufficient for the autotransfusion. If it is necessary to collect more blood, use a new Autotransfusion Unit.
5. Connect a replacement bag for the wound drainage, after 6 hours from the beginning of the surgery.

IMMEDIATE INFUSION OF RECOVERED BLOOD

1. Extract the transfusion set by observing all aseptic procedures and proceed as indicated in the directions for use of the set.
2. The 40 µm filter (10), which is directly inserted in the Autotransfusion Unit, allows eliminating all microaggregates before patient reinfusion.

WARNINGS


- The ATS Bulb Set and its components can be used only if the physician and the assistant staff are aware of the possible implications associated with drainage procedure and autologous blood transfusion.
- It is warmly recommended to observe the "Universal precautions relating to blood and body fluids for treatment of any blood or similar products".
- The staff must always wear gloves.
- Do not handle drainage tubes nor make additional drainage openings.
- If the tube is fixed with suture, it is warmly recommended to avoid possible damages. Do not perform any fixing suture passing through the tube.
- During the drainage tube removal, do not use any forceps or other

Identificazione paziente - Patient identification - Patientenidentifikation - Identification patient - Identificación del paciente - Identificatie patiënt - Identificação do paciente - Pasientidentifikasjon - Patientidentifiaring - Potilastiedot - Patient-id - Αναγνώριση ασθενούς - Identifikační údaje pacienta - Данные пациента

Ora inizio recupero - Recovery start time - Anfangszeit der Erholung - Temps de récupération - Tiempo inicio recuperación - Herstel start tijd - Início da recuperação - Klokkeslett for start av oppsamling - Tidpunkt för start av uppsamling - Talteenoton aloitusaiika - Tidspunkt for start af prøvetagning - Ωρα αρχής ανάκτησης - Godzina rozpoczęcia odzyskiwania - čas zahájení sběru krve - Время начала восстановления

- tools and practise a moderate traction to avoid damages.
- Do not use the drainage tube for more than 29 days.
- When the Autotransfusion Unit is used with a pressure transfusion system (pressure cuff) the staff in charge must be aware of all the warnings and contraindications for this procedure. In any case, make always reference to the warnings and instructions for use of the manufacturer.
- To obtain operational continuity, the ATS Bulb Set must be placed below the level of the patient and with the clamps (13) and the one of the Autotransfusion Unit (14) open.
- The 1,000 ml Autotransfusion Unit (11) and the transfusion set must not be reused. Usually, due to the high capacity of the Autotransfusion Unit, it is not necessary to use a second unit; Redax supplies a package ready for use when required.
- If, by mistake, the closed system is interrupted, do not reinfuse the contaminated blood.
- If an anticoagulant agent is used, it can be infused through the injection point (2).

STERILE - Ethylene oxide sterilized

The product is sterile if the package is intact. Discharge after each single use, do not re-use. Avoid exposure to high temperatures and ultraviolet rays during storage. For discharging and disposal of the device, adopt suitable precautions and act in accordance with the provisions of the law in force for biologically hazardous waste.

MATERIALS USED

SILICONE, PP, PEHD, ABS, PVC, STEEL.
LATEX-FREE MEDICAL DEVICE.

Date of issue of the last version :
see last page : (REV.: XX-XXXX)

LEGEND

- 01 Patient drainage tube
- 02 Injection point
- 03 Red clamp
- 04 Unidirectional valve
- 05 No return valve
- 06 400 cc silicone bulb
- 07 Autotransfusion Unit connection
- 08 Macroaggregate filter
- 09 Connector with pierceable membrane
- 10 40 µm microembolus filter
- 11 Autotransfusion Unit (bag)
- 12 Inox spring for medium negative pressure
- 13 Output clamp
- 14 Autotransfusion Unit clamp
- 15 Cap
- 16 Pierceable cap for diagnostics

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