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## INSTRUCTIONS FOR USE

### DESCRIPTION

The Drentech® Anemos system is a drainage system by gravity which allows, thanks to the concomitant action of two water valves, obtaining effective and balanced drainage in pneumonectomy operations. In these cases there are different requirements: maintain the cavity,

EN

created by removal of the lung, drained and clean; maintain balanced the area of the mediastinum; maintain a slight average vacuum in the thoracic cavity which compensates the tendency of the residual lung to retract. To achieve these objectives, the system is composed of the following parts:

1. A collection chamber (1) made up of two graduated segments. The collection chamber has a maximum capacity of 2,000 cc.
2. A water control chamber for negative pressure (2): this valve allows entry of filtered air into the pleural cavity when the negative pressure inside it exceeds the maximum value of -13 cmH<sub>2</sub>O. The valve, when filled with water, turns red.
3. A water control chamber for positive pressure (3): this valve allows maintaining the maximum level of positive pressure that can develop in the pleural cavity regulated at a value of +1 cmH<sub>2</sub>O. The valve, when filled with water, turns blue.

The different colouring of the liquid in the two valves allows a clear distinction, thus highlighting the different functions.

The system is fitted with a double protection made up of two filters (5) which, located on the two accesses to the outside, allow filtering the air coming into and going out of the system, isolating the latter from the outside environment. The system must exclusively be used by gravity.

#### INSTRUCTIONS FOR USE

The system is suitable for: evacuating air and/or fluid from the thoracic cavity in pneumonectomy operations. Its draining action, which occurs solely by gravity, allows re-establishing and maintaining intrathoracic pressure gradients such as to encourage re-balancing of the mediastinum area, preventing localisation of excessive positive and negative pressure values. The system hence allows obtaining balanced drainage which facilitates postoperative progress of pneumonectomized patients.

**Note:** This device can only be used in safety conditions in the field of applications and using the procedures indicated on this operating instruction sheet for the specific type of product. The manufacturer declines all liability relating to improper use or use different to that indicated.

#### WARNINGS / PRECAUTIONS

- The system must be maintained in vertical position, also during transport, and always below the fluid level of the pleural cavity (at least 50 cm).
- Avoid loops or squashing the patient drainage tube (6).
- Do not obstruct or cover the air inlet and outlet connectors (5).
- Thanks to its particular design, the system is stable and does not require an I.V. pole for use on the floor. The system has a handle (7) for easier transport of the drainage unit. The handle also allows hooking the system to the bed of the patient.
- The system operates solely by gravity. Do not connect the system to a suction source.

#### GENERAL PRECAUTIONS

Do not use if the package has been opened or damaged. Single-use product. Reuse may lead to alteration of performance and risks of cross-contamination.

#### PREPARING THE SYSTEM

The system is packed in a single wrapping with an internal sheet. To use the product, grip the lower end of the package and pull in the direction of the opening. Extract the product and unfold the internal sheet which can be laid on the floor before positioning the system on it.

##### a. *Filling the positive pressure control chamber (3).*

1. Prepare a syringe.
2. Remove the red filter cap from the connector (4) before filling the control chamber (3).
3. Insert the syringe in the connector (4) and fill the control chamber (3) with sterile water up to the broken line, corresponding to about 60 cc. The water will turn blue thanks to a dye in the chamber.
4. Remove the syringe and place the red cap on the connector.
5. During operation the oscillations synchronous with the breaths of the patient can be observed. These oscillations reflect the variations in the endopleural pressure.
6. The bubbles appearing in the control chamber show the amount of air loss of the patient.

##### b. *Filling the negative pressure control chamber (2).*

1. Remove the closing cap (8) of the filling connector.
2. Fill the control chamber (2) with sterile water up to the broken line, corresponding to about 120 cc. The water will turn red thanks to a dye in the chamber.
3. Reposition the closing cap.
4. The graduated scale indicates the value 13 cmH<sub>2</sub>O, corresponding to the negative pressure level beyond which air input from the valve commences. When this event occurs, bubbles can be seen rising in the widest section of the valve.
5. During operation the oscillations synchronous with the breaths

of the patient can be observed. These oscillations reflect the variations in the endopleural pressure.

#### c. *Connecting the drainage tube*

1. Connect the drainage tube (6) to the thoracic catheter of the patient by means of the connector (6). At this point the system is ready to operate by gravity.

#### d. *Drainage sampling*

In order to evaluate the quality of the drained fluids to ascertain any presence of sepsis, the sample can be taken directly from the patient drainage tube (6). For this purpose, temporarily shape the tube, creating a loop, and obliquely insert the needle of a syringe.

**Note:** Do not perforate the patient tube with a needle larger than 18 Gauge (1.24 mm).

#### e. *System replacement*

To replace the system when it has reached its maximum collection capacity, proceed as follows:

1. Prepare a new system as described in the previous points.
2. Close the drainage upstream of the connector (9) with a clamp.
3. Disconnect the system in use and connect the new one.
4. Open the clamp to restore drainage.

#### CONTRAINDICATIONS

There are no known contraindications.

#### STERILE - Ethylene oxide gas sterilized.

The product is sterile unless the package has been opened or damaged. Discard after each single use, do not re-use. Avoid exposure to high temperatures and ultraviolet rays during storage. For discarding and disposal of the device, adopt suitable precautions and act in accordance with the provisions of the law in force for biologically hazardous waste.

**Note:** The liquid contained in the collection chambers of the drainage system must absolutely not be used for infusion.

#### MATERIALS USED

PS, PP, S-E/B-S.  
LATEX-FREE MEDICAL DEVICE.

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