REDAX®

# DRENTECH™ COMPACT DRENTECH™ DOUBLE

IT

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### DESCRIZIONE

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## EN

#### DESCRIPTION

The systems of the Drehtech™ Compact family are single-use thoracic drainage devices for adult and paediatric use equipped with a water seal valve for suction control. They can be used by suction or gravity and are available in single or double tube version.

The drainage systems are composed of the following parts:


1. A collection chamber (1) divided into three graduated segments with a sampling port on the rear of the drain (10). The collection chamber has a maximum capacity of 2100 ml. In the version with double tube, the collection chamber is divided into two independently graduated segments to allow separation of the drainage fluids.
2. A water seal valve (2) which allows air evacuation from the thoracic cavity of the patient preventing backflow. The seal valve (2) has a graduated scale.
3. Two water columns (3) to control suction, which allows regulating the negative pressure by means of the height of the water level inside them (from -5 cmH<sub>2</sub>O to -25 cmH<sub>2</sub>O).

#### INDICATIONS FOR USE

The drainage systems are indicated for:

- Evacuating air and/or fluids from the thoracic cavity or the mediastinum preventing their build-up.
- Re-establishing and maintaining normal pressure gradients in the thoracic cavity.
- Facilitating complete pulmonary expansion to restore normal respiratory dynamics.

**Note:** This device can be used in safe conditions only in the field of application and according to the procedures indicated in

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se siste side : (REV.: XX-XXXX)

this instruction leaflet in correspondence to the type of product. The manufacturer declines all responsibility for improper use or different from that indicated.

#### GENERAL WARNINGS

- The product is sterile if the package is intact. Do not use if the package has been opened or damaged.
- Dispose of the device after each single use; do not reuse.
- Avoid exposure to high temperatures and ultraviolet rays during storage.

#### WARNINGS/PRECAUTIONS

- **The drainage system must be held in vertical position and always below the fluid level of the pleural cavity (at least 50 cm).**
- **Avoid looping or squashing the patient drain tube (9).**
- **Do not obstruct or cover the positive pressure relief valve (5).**
- **Do not connect the system to the patient with the suction source already activated.**
- **Thanks to its particular design, the system is stable and does not require a tube support for use on the floor. The system has a handle (12) to easier the transport of the drainage unit. The handle also allows hooking the system to the bed of the patient.**
- **Do not clamp the connection tube to the patient for an extended period of time during drainage and/or transport.**
- **The fluid contained in the collection chambers of the drainage system may absolutely not be reinfused, but must be disposed of with the drainage system in accordance with the regulations in force on biologically hazardous waste.**

#### PREPARING THE DRAINAGE SYSTEM

The system is packed in double sterile wrapping. Remove the wrapping using an aseptic technique.

##### a. Filling the water seal valve

1. Fill a syringe with sterile water or saline solution. The needle does not have to be used.
2. Insert the syringe in the connector (7) and fill the seal valve (2) up to the broken line (about 45ml as also indicated on the front of the device). The water will turn blue thanks to a dye contained in the valve.
3. Remove the syringe: the connector (7) is equipped with a valve that automatically closes its access when the syringe is removed.
4. At this point, the seal valve is ready.

##### b. Filling the suction regulating columns

1. Remove the protective cap (4) of the control columns (3).
2. Fill the control columns (3) with sterile water or saline solution up to the desired suction level (max. 25 cmH<sub>2</sub>O). The water contained in the chamber will turn blue. The amount of liquid required depends on the desired pressure value; the approximate values corresponding to the pressures indicated are given below.

H<sub>2</sub>O volume suction level:

10 cmH<sub>2</sub>O 160 ml

15 cmH<sub>2</sub>O 250 ml

20 cmH<sub>2</sub>O 320 ml

25 cmH<sub>2</sub>O 390 ml

3. Wait until the water level in both control columns (3) has aligned and corresponds to the desired suction level.
4. Close the cap (4) and close the front valve (11) as shown in the figure.
5. The water level in the control columns (3) determines the negative pressure values applied to the patient. The correct suction level must be chosen by a physician.

**WARNING:** it is important to check that the front valve (11) is correctly closed before connecting the suction source and starting to use the device. If the valve is not correctly closed, the columns would be unbalanced until partial fluid suction to the centralized source. To restore the correct levels, see the "Troubleshooting" table.

##### c. Connecting the drain tube

Connect the drain tube (9) to the thoracic catheter of the patient and then connect the other end to the collection chamber by means of the lock-in connector (13). In the version with four chambers, there are two connection tubes for connection to the thoracic catheters, identifiable by the numbering printed directly on the tube.

If you intend to use the system to treat neonatal patients, use the dedicated tubing for connection to the built-in Heimlich valve (code Redax 10564).

#### DEVICE OPERATION

##### a. Operation by suction

1. Connect a connection tube to the connector (6) to connect the drainage system to the suction source.

**Note:** The drainage system is equipped with an automatic "self-regulating" device (16) which allows obtaining moderate gurgling in the pressure regulating columns (3) in any vacuum condition available in hospital departments and without needing any intervention by the staff. For this reason, the vacuum can directly be connected to the system without interposing other valves and/or regulators.

##### WARNINGS:

- Periodically check that suction in the control chamber (3) is continuous with moderate gurgling.
- Suction must be suspended in the following cases: while restoring the water levels or when you want to replace the system. After each suspension, check that suction is restored.
- Evaporation, especially at high temperatures, may cause the liquid level in the control columns (3) and in the water seal valve (2) to drop. This evaporation is normally negligible thanks to the "self-regulating" device. It is in any case advisable to periodically check the liquid level and top it up if necessary.

##### b. Operation by gravity

The drainage system can also be used by gravity by detaching the suction tube connected to the vacuum source from the connector (6). The automatic positive pressure relief valve (5) allows operation by gravity without requiring further operations.

##### c. Intrathoracic pressure diagnostics

During patient drainage, there may be variations in intrathoracic pressure, qualitatively viewable as oscillations of the liquid in the graduated column of the seal valve (2). The actual pressure applied to the patient is obtained from the sum of the pressure determined by the regulating columns (3) and from that read on the column of the seal valve (2). In case of operation by gravity, the first factor is null.

The oscillations synchronous with the patient's respiratory acts reflect the intrathoracic pressure changes and may be an indication of complete lung expansion or drain tube occlusion. In case of doubt, the physician has to do the necessary patency checks.

##### d. Patient air leakage diagnostics

Air bubbles in the lower chamber of the seal valve (2) are an indication of air leakage and, consequently, incomplete lung expansion. A reduction in the amount of bubbles over time indicates a gradual reduction of the leakage and better lung expansion. When lung expansion is complete, there will be no bubbles in the seal valve.

**WARNING:** If the phenomenon persists, ask a physician to intervene after carefully checking the connections.

##### e. High negative pressure relief valve

Excessive negative pressure (e.g. caused by the patient's respiratory dynamics or by the effect of patient tube "milking") can be seen from the height of the water column in the seal valve (2).

Generally, this negative pressure spontaneously reduces following patient air and fluid leakage. If you want to manually reduce this pressure, insert the special pin (8) in the connector (7) equipped with a dedicated valve; the latter opens allowing a modest amount of air to flow into the collection chamber thus reducing the negative pressure inside it. This reduction can be controlled based on the height of the water in the column of the seal valve. Release the pin when the desired negative pressure level has been reached and in any case not beyond level 0 (zero); the valve in the connector (7) will automatically close.

In cases of high negative pressure, the float valve protects the water seal valve in the presence of high negative pressure preventing liquid backflow. The water induces closure of the float valve when the negative pressure becomes excessive; the valve reopens when the negative pressure is reduced.

**WARNING:** do not press the high negative pressure relief valve when the system is operating by gravity or when suction is not active.

##### f. Reducing the suction level

The suction level can be reduced at any time if deemed necessary. The device allows easily and quickly reducing the suction level to half its original value. To carry out this operation, just open the protective cap (4). You will immediately notice that the characteristic bubbling disappears in the left-hand column while it remains

unchanged in the right-hand column. That way, one of the two regulating columns is excluded thus reducing suction to about half the value originally set.

To restore the initial suction, just reposition the cap (4). The operations described can be carried out while the device is in normal operation and suction does not have to be suspended. Evaporation, especially at high temperatures, or considerable air leakage from the patient may cause the liquid level in the seal valve to drop. It is advisable to periodically check the liquid level and top it up if necessary.

#### g. Transporting the device

The device should preferably be transported in vertical position and below the level of the patient.

If necessary, the device can be placed in horizontal position as shown in Fig. 1. In this position, avoid mixing the liquids between the various chambers.

#### IMPORTANT

- In horizontal position, the water seal is compromised and the patient tube must therefore be clamped. This operation must be assessed by the responsible physician in order to evaluate the risks that may arise from clamping the patient.
- To avoid prolonged suspension of the drain with the risk of hypertensive pneumothorax, the clamp must be re-opened as soon as transport is complete. In any event, the clamp must not remain closed for more than 15 consecutive minutes.

#### h. Replacing the unit

Replace the drainage system when the capacity of the collection chamber has been used up. To carry out this operation:

1. Prepare the new unit following the steps above so that it is ready for use.
2. Suspend suction and close the patient tube by means of the sliding clamp provided.
3. Disconnect the patient tube from the unit in use by means of the lock-in connector and replace it with the new one.
4. Re-open the clamp on the patient tube and restore suction.

#### DISPOSAL

Adopt adequate precautions for elimination and disposal of the device and comply with the provisions of the laws in force on biologically hazardous waste.

#### CONTRAINDICATIONS

There are no known contraindications.

#### SINGLE-USE DEVICE.

Reuse may lead to alterations in performance and cross-contamination risks.

#### STERILE - Sterilized by ethylene oxide

#### MATERIALS USED

Polystyrene (PS), Polypropylene (PP), Elastomer (S-E/B-S), Polyethylene (PEHD), ABS.

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