cod. 10564

Redax S.p.A.  
Via Galileo Galilei, 18  
Poggio Rusco (MN) Italy

COD. 00204

REV.: 05-2021

IT

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DESCRIZIONE

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EN

#### DESCRIPTION

The Drentech™ Variant and Drentech™ Dune systems are single-use thoracic drainage devices for adult and paediatric use equipped with a collection chamber, seal valve and mechanical regulator for suction control. The Variant systems have a water seal valve and the Dune systems a mechanical seal valve. The latter operate completely "dry" and use of water is optional solely for the purpose of diagnosing any patient air leakage. All the devices can be used by suction or gravity and are available in single or double tube version.

The systems are composed of the following parts:

1. A collection chamber (1) divided into three graduated segments with a sampling port on the inlet of the drain (10). The collection chamber has a maximum capacity of 2100 ml. In the version with double tube, the collection chamber is divided into two independently graduated segments to allow separation of the drainage fluids.
2. A seal valve (2). The Variant systems have a water seal valve and the Dune systems a mechanical seal valve.
3. A mechanical suction control regulator (3). The optimal suction flow setting for valve operation can be viewed by means of the dedicated float indicator (4). The suction level can be established by means of a graduated knob (11) in a range between -5 and -40 cmH<sub>2</sub>O.

Fig. 1




Fig. 2

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INDICATIONS FOR USE

The drainage systems are indicated for:

- Evacuating air and/or fluids from the thoracic cavity or the mediastinum preventing their build-up.
  - Re-establishing and maintaining normal pressure gradients in the thoracic cavity.
  - Facilitating complete pulmonary expansion to restore normal respiratory dynamics.
- Note:** This device can be used in safe conditions only in the field of application and according to the procedures indicated in this instruction leaflet in correspondence to the type of product. The manufacturer declines all responsibility for improper use or different from that indicated.

GENERAL WARNINGS

- The product is sterile if the package is intact. Do not use if the package has been opened or damaged.
- Dispose of the device after each single use; do not reuse.
- Avoid exposure to high temperatures and ultraviolet rays during storage.

WARNINGS/PRECAUTIONS

- **The drainage system must be held in vertical position and always below the fluid level of the pleural cavity (at least 50 cm).**
- **In the Variant systems, the water seal is compromised in horizontal position. It is advisable to use a clamp on the patient tube for as long as the system remains in horizontal position.**
- **In the Dune systems, the seal remains active also in horizontal position so that it is not necessary to clamp the patient tube.**
- **Avoid looping or squashing the patient drain tube (9).**
- **Do not obstruct or cover the positive pressure relief valve (5).**
- **Do not connect the system to the patient with the suction source already activated.**
- **Do not clamp the connection tube to the patient for an extended period of time during drainage and/or transport.**
- **The fluid contained in the collection chambers of the drainage system may absolutely not be reinfused, but must be disposed of with the drainage system in accordance with the regulations in force on biologically hazardous waste.**

PREPARING THE DRAINAGE SYSTEM

The systems are supplied in double sterile packaging. Unpack the product using an aseptic technique.

a. Filling the water seal valve

- *Only for VARIANT systems*

**NOTE:** The Variant drainage systems have a seal valve (2) which requires filling with liquid in order to function. Only in this way is creation of the water seal that prevents air backflow to the pleural cavity guaranteed.

1. Fill a syringe with sterile water or saline solution. The needle does not have to be used.
  2. Insert the syringe in the connector (7) and fill the seal valve (2) up to the broken line (about 45ml as also indicated on the front of the device). The water will turn blue thanks to a dye contained in the valve.
  3. Remove the syringe: the connector (7) is equipped with a valve that automatically closes its access when the syringe is removed.
  4. At this point, the seal valve is ready.
- *Only for DUNE systems*

In the Dune systems, the seal valve is mechanical and hence does not require filling with water.

This operation can optionally be carried out for air leakage diagnostics (see the paragraph "Air leakage diagnostics" in this instruction leaflet).

b. Connecting the drain tube

Connect the drain tube (9) to the thoracic catheter of the patient and then connect the other end to the collection chamber by means of the lock-in connector (13). In the version with four chambers, there are two connection tubes for connection to the

thoracic catheters, identifiable by the numbering printed directly on the tube.

If you intend to use the system to treat neonatal patients, use the dedicated tubing for connection to the built-in Heimlich valve (code Redax 10564).

DEVICE OPERATION

a. Operation by suction

Connect a connection tube to the connector (6) to connect the drainage system to the suction source.

b. Regulating the suction flow

The drainage system is equipped with an automatic "self-regulating" device (14) which allows obtaining optimal functioning of the regulating valve in any vacuum condition available in hospital departments and without needing any intervention by the staff. For this reason, the vacuum can directly be connected to the system without interposing other valves and/or regulators. After connecting the suction source, the red float indicator (4) will be positioned and clearly visible in the area identified by the transparent window on the front (see Figure 2).

WARNINGS:

- If the float indicator is not visible as it is in the area below the window, the regulator may be unable to reach the higher suction values. Check the centralized suction source in order to obtain optimal positioning of the float indicator.
- **The proper working of the mechanical valve is guaranteed, whenever the red floater is visible through its designated window.**

c. Regulating the suction

- The suction applied to the patient can be regulated from -5 cmH<sub>2</sub>O to -40 cmH<sub>2</sub>O (1 cmH<sub>2</sub>O = -0.1 kPa = -0.1 mbar = -0.014 Psi).
- To determine the suction level, turn the regulating knob (11) until reaching the desired value.
- The negative pressure immediately stabilizes and the system is ready for use.

d. Operation by gravity

The drainage system can also be used by gravity by detaching the suction tube connected to the vacuum source from the connector (6). The automatic positive pressure relief valve (5) allows operation by gravity without requiring further operations.

e. Intrathoracic pressure diagnostics (Variant systems only)

During patient drainage, there may be variations in intrathoracic pressure, qualitatively viewable as oscillations of the liquid in the graduated column of the seal valve (2). The actual pressure applied to the patient is obtained from the sum of the pressure determined by the mechanical regulating valve (3) and from that read on the column of the seal valve (2). In case of operation by gravity, the first factor is null.

The oscillations synchronous with the patient's respiratory acts reflect the intrathoracic pressure changes and may be an indication of complete lung expansion or drain tube occlusion. In case of doubt, the physician has to do the necessary patency checks.

**NOTE:** in the Dune systems, the presence of the mechanical seal valve does not allow viewing the oscillations described above and therefore neither qualitative evaluation of the intrathoracic pressure.

f. Patient air leakage diagnostics

This function is available in both versions of the device after filling the lower chamber of the seal valve (2). Air bubbles in this chamber are an indication of air leakage and, consequently, incomplete lung expansion. A reduction in the amount of bubbles over time indicates a gradual reduction of the leakage and better lung expansion. When lung expansion is complete, there will be no bubbles in the seal valve.

**WARNING:** If the phenomenon persists, ask a physician to intervene after carefully checking the connections. Evaporation, especially at high temperatures, or considerable air leakage from the patient may cause the liquid level in the seal valve to drop. It is advisable to periodically check the liquid level and top it up if necessary.

g. High negative pressure relief valve

Excessive negative pressure (e.g. caused by the patient's respiratory dynamics or by the effect of patient tube "milking") can be seen from the height of the water column in the seal valve (2).

Generally, this negative pressure spontaneously reduces following patient air and fluid leakage. If you want to manually

reduce this pressure, insert the special pin (8) in the connector (7) equipped with a dedicated valve; the latter opens allowing a modest amount of air to flow into the collection chamber thus reducing the negative pressure inside it. This reduction can be controlled based on the height of the water in the column of the seal valve. Release the pin when the desired negative pressure level has been reached and in any case not beyond level 0 (zero); the valve in the connector (7) will automatically close.

In cases of high negative pressure, the float valve (4) protects the water seal valve in the presence of high negative pressure preventing liquid backflow. The water induces closure of the float valve when the negative pressure becomes excessive; the valve reopens when the negative pressure is reduced.

**WARNING:** do not press the high negative pressure relief valve when the system is operating by gravity or when suction is not active.

#### h. Transporting the device

Thanks to its particular design, the system is stable and does not require an IV pole to use it on the floor. The system has a handle (12) to facilitate transport of the drainage unit. The handle moreover allows hooking the system onto the patient's bed.

The device should preferably be transported in vertical position and below the level of the patient.

If necessary, the device can be placed in horizontal position as shown in Fig. 1. In this position, avoid mixing the liquids between the various chambers.

#### IMPORTANT

##### For VARIANT systems

- In horizontal position, the water seal is compromised and the patient tube must therefore be clamped. This operation must be assessed by the responsible physician in order to evaluate the risks that may arise from clamping the patient.
- To avoid prolonged suspension of the drain with the risk of hypertensive pneumothorax, the clamp must be re-opened as soon as transport is complete. In any event, the clamp must not remain closed for more than 15 consecutive minutes.

#### i. Replacing the unit

Replace the drainage system when the capacity of the collection chamber has been used up. To carry out this operation:

1. Prepare the new unit following the steps above so that it is ready for use.
2. Suspend suction and close the patient tube by means of the sliding clamp provided.
3. Disconnect the patient tube from the unit in use by means of the lock-in connector and replace it with the new one.
4. Re-open the clamp on the patient tube and restore suction.

#### DISPOSAL

Adopt adequate precautions for elimination and disposal of the device and comply with the provisions of the laws in force on biologically hazardous waste.

#### CONTRAINDICATIONS

There are no known contraindications.

#### SINGLE-USE DEVICE.

Reuse may lead to alterations in performance and cross-contamination risks.

#### STERILE - Sterilized by ethylene oxide

#### MATERIALS USED

Polystyrene (PS), Polypropylene (PP), Elastomer (S-E/B-S), Polyethylene (PEHD), ABS.

#### MEDICAL DEVICE NOT MADE FROM NATURAL LATEX PVC-FREE PACKAGING

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