Quantità di Sangue - Blood quantity - Blutmenge - Quantité de sang - Cantidad de sangre - Hoeveelheid bloed - Quantidade de sangue - Blodmengde - Blodmängd - Verimäärä - Blodmængde - Ποσότητα αίματος - Количество крови





Annotazioni - Remarks - Bemerkungen - Remarque - Observaciones - Opmerkingen - Anotações - Merknader - Anmärkningar - Huomautuksia - Bemærkninger - Σημειώσεις - Замечания

Identificazione paziente - Patient identification - Patientenidentifikation - Identification patient - Identificación del paciente - Identificatie patiënt - Identificação do paciente - Pasientidentifikasjon - Patientidentifisering - Potilastiedot - Patient-id - Αναγνώριση ασθενούς - Данные пациента

Ora inizio recupero - Recovery start time - Anfangszeit der Erholung - Temps de récupération - Tiempo inicio recuperación - Herstel start tijd - Início da recuperação - Klokkeslett for start av oppsamlng - Tidpunkt for start av uppsamlng - Talteenoton aloitusaika - Tidspunkt for start af prøvetagning - αρχής ανάκτησης - Время начала восстановления

COD. 00181

REV.: 05-2017

Redax S.p.A.
Via Galileo Galilei, 18
Poggio Rusco (MN) Italy

REDAX™

IT

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DESCRIZIONE

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EN

DESCRIPTION

The Drentech™ Emotrans Dune systems are single-use thoracic drainage devices fitted with collection chambers, seal valve and mechanical regulator for suction control. The seal valve is made by means of a silicone membrane designed so that the valve is normally closed. In this way, the seal is always assured in all situations. The Dune systems operate in a completely "dry" manner and the use of water is optional with the only purpose to diagnose possible patient air loss. Drainage systems can be used both by suction or gravity and are available in both adult and paediatric versions, with three chambers (single tube). The systems are composed of the following parts:

1. A collection chamber (1) with a maximum capacity of 2.150 ml. The first section features a filter, with 200 µm porosity and large surface for the preliminary filtering of blood and macroaggregates elimination.
2. A seal valve (2). This valve is mechanical, made by means of a silicone membrane and seal guaranteed when it is closed. Therefore, the "seal" is always guaranteed in any condition and there is no risk of blood contamination. The water in the chamber (2) is used exclusively to detect air loss and does not constitute any underwater seal.
3. A mechanical-type suction control regulator (3). The setting of the optimal suction flow for valve operation can be viewed by means of the floating indicator (4). The suction level can be fixed by means of a graduated knob (12), in an interval between -5 and -40 cmH₂O.
4. Back connector (8) for the connection to the Redax blood recovery bag, fitted with a connector with a pierceable membrane.
5. A label, situated on the side of the system, on which it is possible to put patient data and the starting time of autologous blood recovery.

INSTRUCTIONS FOR USE

The Emotrans systems are intended for:

- Evacuating air and/or fluids from the chest cavity or the mediastinum, preventing their accumulation;
- Re-establishing and maintaining normal pressure gradients in the thoracic cavity;
- Facilitating complete pulmonary expansion to restore normal respiratory dynamics;
- Collecting blood during postoperative phase through a preliminary filtering to eliminate microaggregates and surgical residuals.

The system is set for collected blood recovery, for immediate autotransfusion, thanks to the possibility of connection with the Redax blood recovery bag. Moreover, the system can be connected, by means of a suitable optional connection tube, to a separating device for the possible concentration and/or lavage of the erythrocytes before the reinfusion.

GENERAL WARNINGS

- Do not use if the package has been opened or damaged.
- This device can only be used in safe conditions in the field of applications and using the procedures indicated on this operating instruction sheet for the specific type of product. The manufacturer declines all liability relating to improper use or use different from that indicated.
- This device must be used by health care personnel trained for the use of thoracic drainage systems and blood recovery systems, with the supervision of a physician.

WARNINGS/PRECAUTIONS

- The Emotrans system must be maintained in a vertical position, also during transport, and always below the fluid level of the pleural cavity (at least 50 cm). The system can be put in a horizontal position, as shown in Figure 1, for a short time (ex. during patient transfer), without fluid blending between the different chambers or fluid release.
- The seal remains active also in horizontal position so it is not necessary to clamp the patient tube.
- In the event that the system overturns accidentally, the mechanical seal valve does not let the fluids container in the system blend with the collected blood. This latter remains unchanged and it is available for a subsequent reinfusion.
- Avoid kinks or squashing the patient drainage tube (9).
- Do not obstruct or cover the positive pressure release valve (5).

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- Do not connect the system to the patient when the suction is already into operation.
- Do not clamp the patient connection tube for a long time during the drainage and/or the transport.
- To use the system to collect and recover autologous blood, together with the Redax blood recovery bag, refer to the instructions for use of the latter.
- To use the Redax optional connection tube, to connect to a separating device, refer to the corresponding paragraph of this manual of instructions for use.
- As recovered blood is kept at room temperature, the time between the starting of the collection and the starting of the reinfusion must be within and not later than 6 hours.

PREPARING THE SYSTEM

The systems are supplied in double sterile packaging. Unpack the product using an aseptic technique.

a. Preparing the seal valve

In the Emotrans systems, the seal valve is made by means of a silicone membrane so it does not require to be filled with water. This operation can be carried out, in an optional manner, to diagnose air loss (see the relevant paragraph in this instructions sheet). In this way, the valve guarantees air flow unidirectionality and the protection of the blood collected in the chamber in all situations.

b. Connecting the drain tube

Connect the drain tube (9) to the thoracic catheter of the patient and then connect the other end to the collection chamber by means of the lock-in connector (11).

DEVICE OPERATION

a. Operation by suction

Connect a connection tube to the connector (6) to connect the drainage system to the suction source.

b. Regulating the suction flow

The drainage system is equipped with an automatic "self-regulating" device (14) which allows obtaining optimal functioning of the regulating valve in any vacuum condition available in hospital departments and without needing any intervention by the staff. For this reason, the vacuum can directly be connected to the system without interposing other valves and/or regulators. After connecting the suction source, the red float indicator (4) will be positioned and clearly visible in the area identified by the transparent window on the front (see Figure 2). This position identifies the optimal suction flow.

WARNINGS:

- If the float indicator is not visible as it is in the area below the window, the regulator may be unable to reach the higher suction values. Check the centralized suction source in order to obtain optimal positioning of the float indicator.
- **The proper working of the mechanical valve is guaranteed, whenever the red floater is visible through its designated window.**

c. Regulating the suction

- The suction applied to the patient can be regulated from -5 cmH₂O to -40 cmH₂O (1 cmH₂O ≈ -0.1 kPa ≈ -0.1 mbar ≈ -0.014 Psi).
- To determine the suction level, turn the regulating knob (12) until reaching the desired value.
- The negative pressure immediately stabilizes and the system is ready for use.

d. Operation by gravity

The drainage system can also be used by gravity by detaching the suction tube connected to the vacuum source from the connector (6). The automatic positive pressure relief valve (5) allows operation by gravity without requiring further operations.

e. Diagnostics of patient air loss

In the Emotrans Dune systems, the observation of air loss can be carried out through the lower chamber filling (2). As the seal is guaranteed by the silicone membrane mechanical valve this filling is optional and depends only on the physician need to evaluate the above-mentioned air loss. To fill the lower chamber, act as follow:

1. Fill a syringe with sterile water or physiological solution. It is not necessary to use the needle.
2. Insert the syringe in the connector (7) and fill the chamber (2), up to the broken line (about 45 ml). The water will turn blue thanks to a dye present in the valve.
3. Remove the syringe: the connector (7) is provided with a valve that stops access automatically when the syringe is removed.

WARNINGS

- The presence of air bubbles in the lower chamber of the seal valve (2) indicates an incomplete lung expansion. The reduction of the quantity of these bubbles in the course of time indicates a progressive lung expansion. When this expansion is complete, it will be possible to verify the absence of bubbles

in the seal valve.

- Extensive bubbling in the seal valve (2) can be indicative of a serious pneumothorax or air loss in the connections. Carry out a careful checking and, if the phenomenon persists, call in a physician.

f. Reducing the suction level

The suction level can be reduced at any time if necessary on specification of the physician in charge. To carry out this operation, act on the regulating knob (12), until reaching the desired suction pressure. It is not necessary to carry out any suction flow regulation.

g. Transporting the device

Thanks to its particular design, the system is stable and does not require an IV pole to use it on the floor. The system has a handle (10) to facilitate transport of the drainage unit. The handle moreover allows hooking the system onto the patient's bed.

The device should preferably be transported in vertical position and below the level of the patient.

If necessary, the device can be placed in horizontal position as shown in Fig. 1. In this position, avoid mixing the liquids between the various chambers.

h. Replacing the unit

Replace the drainage system when the capacity of the collection chamber has been used up. To carry out this operation:

1. Prepare the new unit following the steps above so that it is ready for use.
2. Suspend suction and close the patient tube by means of the sliding clamp provided.
3. Disconnect the patient tube from the unit in use by means of the lock-in connector and replace it with the new one.
4. Re-open the clamp on the patient tube and restore suction.

i. Collection and recovery of autologous blood

For bag connection operations and bag use, refer to the respective operating instructions.

WARNINGS

- As the recovered blood is kept at room temperature, the time intervening between the beginning of the collection and the beginning of the reinfusion must not exceed 6 hours: refer to the details shown on the label. In the event that this information is not reported, refer to the time of the end of the surgery.
- Normally, it is not necessary to add coagulating substances. If the physician decides to use anticoagulant, on the basis of considerations concerning the clinical picture of the patient, it is advisable to use ACD-A (in a dosage of 40-60 ml for each unit of collected blood or in an equivalent ratio of 1:7) because it is better metabolised. In patients with preserved hepatic functionality, the ACD-A anticoagulant solution does not require specific antagonists if the dosage of 150-180 ml is not exceeded. The calcium solution administration antagonises the effects of a possible excessive anticoagulant administration.

The presence of a tank under the infusion point lets a part of the infused anticoagulant to remain contained. In this way, the drained blood flows directly in this tank blending immediately with the anticoagulant that holds up the coagulation stimulus.

CONTRAINDICATIONS

No known contraindications exist for drainage use. For the recovery and reinfusion of autologous blood, see the Redax blood recovery bag operating instructions. In general, it is best to follow the "Universal precautions relating to blood and body fluids".

DISPOSAL

Adopt adequate precautions for elimination and disposal of the device and comply with the provisions of the laws in force on biologically hazardous waste.

SINGLE-USE DEVICE.

Reuse may lead to alterations in performance and cross-contamination risks.

STERILE - Ethylene oxide sterilized.

MATERIALS USED

Polystyrene (PS), Polypropylene (PP), Elastomer (S-E/B-S), Polyethylene (PEHD), ABS, PVC, LATEX-FREE MEDICAL DEVICE. PVC-FREE PACKAGE.

Date of issue of the last version :

see last page : (REV.: XX-XXXX)