
Fig. 13

REDAX™

IT

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EN

**DECLARATION OF LIABILITY**

This user manual is supplied in printed form and accompanies the device.

This document is the exclusive property of Redax S.p.A. all rights reserved. Any disclosure, reproduction or transfer of content to third parties is prohibited without the prior authorization of the Company.

Redax S.p.A. denies all liability with respect to damage to people or property due to improper use of this product and failure to observe the indications, warnings, instructions and precautions in this user manual.

In the event of a serious accident related to the medical device, the user must notify the manufacturer by e-mail [vigilance@redax.it](mailto:vigilance@redax.it) and the competent authority of the state where the event occurred.

**MANUFACTURER**

**Redax S.p.A.**  
Via Galileo Galilei, 18, 46025 Poggio Rusco (MN) Italia  
P.IVA: 02556750368  
[info@redax.it](mailto:info@redax.it)  
[www.redax.it](http://www.redax.it)

**1. DEVICE DESCRIPTION**

The 'UNICO' device is a percutaneous drain set available in different versions and in different diameters. The set consists essentially of the following components:

1. A main body with two entrances in the Forty and Base (Fig.1a) versions or three entrances (Fig.1b) in all other versions.
2. The first access point (1), located posteriorly, is for the introduction of the Verres needle or the atraumatic mandrel where provided for by the positioning procedure. The two remaining access points (2) are equipped with luer-lock connections and allow the connection of drain devices such as syringes, collection bags, suction systems, etc. in the manner described below in this instruction sheet. On the main body there is a selector switch (3) which allows, by means of a 90° rotation, either gravity drainage or the injection of liquids such as washing solutions. In the versions 'without the extension line', the drain catheter is permanently attached to the main body (4). In the 'with extension' versions, the catheter is connected to the main body by means of a butterfly connection (5) that allows the interposition of an extension line (6). In all versions, except Unico Base, the catheter has holes and

depth markings.

- All versions feature a Verres-type needle (7) for introducing the catheter. This needle has a mandrel with an atraumatic tip that protrudes from the end of the needle. This mandrel can slide longitudinally until the cutting profile of the needle itself is uncovered during insertion, and then return to its extended position by means of a spring. The needle is equipped with a handle (8) in which there is a visual indicator that informs the user of the position of the protective atraumatic mandrel. The Verres needle can be introduced through the first access (1) or, for the Forty version, through the special hole in the catheter (Fig.2).
- The devices are available in a 'standard' version in which the catheter is straight and in a 'multi' version in which the catheter has a distal 'pig-tail' end (Fig.3).

#### STERILE - STERILISED WITH ETHYLENE OXIDE.

The product is sterile if the packaging is intact.

The device is disposable and therefore requires no cleaning or maintenance.

#### STORING

The device must be stored according to the label: temperature limit (Min/Max) 0°C/60°C.

## 2. DIRECTIONS FOR USE

The 'UNICO' system is primarily intended to create a minimally invasive access to the patient's thoracic and abdominal cavities in order to evacuate liquids and air accumulated there, or a combination of the two. The device can be used for gravity or suction drainage, connected to collection systems. The device can be used for the infusion of liquids such as physiological flushing solutions or talcum powder nebulised in sterile injectable solutions.

#### PATIENT GROUP

The device is indicated for use on adult patients.

#### USER GROUP

All products are intended for use by qualified medical personnel, such as surgeons and nurses.

#### MEDICAL CONDITIONS

The fields of application for UNICO drains are: treatment of pneumothorax (primary and secondary), malignant or chronic pleural effusions, empyema and abdominal ascites.

#### CONTRAINDICATIONS

Coagulopathies and other blood diseases.

#### INSTRUCTIONS

- Do not re-sterilise the device.
- Disposable product. Reuse can lead to altered performance and risk of cross-contamination.
- Do not use the device beyond the expiry date indicated on the packaging.
- Avoid exposure to high temperatures and direct UV light during storage.
- No modification of the device and/or its parts is permitted.
- This device may only be used safely in the field of application and in the manner specified in this instruction sheet, corresponding to the product type. The manufacturer accepts no liability for improper use or use other than that indicated.
- Check that the visual indicator of the mandrel position runs freely and that the corresponding indicator is green when the mandrel is in the extended position. If, during initial compression, the green indicator on the outer surface of the chest wall does not disappear, stop the needle introduction procedure and replace the device.
- The device is designed to be used with Redax products and accessories. If different products and accessories are used, it is the doctor's responsibility to check their compatibility before use.
- Do not use the device for more than 30 days.
- The duration of the drainage procedure must be determined by the doctor based on clinical considerations.

- Do not use the device if the casing or any of its parts are damaged.
- Do not manipulate the drainage in any way or drill additional holes using cutting equipment.
- In the case of suture fixation, carefully avoid suturing the tube itself or cutting it with needles or sharp objects. Avoid excessively tight sutures as this may cause drain rupture.
- After positioning the drain, ensure that the perforated section is completely included within the wound and that the connections are tight.
- Drain should be removed manually by gently pulling, avoiding sudden manoeuvres. Avoid using metal instruments such as pliers or devices that could cause drain rupture.
- After removal, check that the drainage set is complete as this is essential in order to exclude the presence in the cavity of drainage fragments caused by accidental tearing.
- Improper use could expose the patient and the operator to various dangers.
- Keeping the drain in situ for very long periods, could make its removal difficult. Care should be taken during removal operations.

## 3. CLINICAL BENEFIT

The clinical benefit associated with UNICO percutaneous drains is to create a minimally invasive access to the patient's thoracic or abdominal cavities in order to evacuate accumulated fluids (air and pathological fluids) and alleviate related symptoms. The main feature of the device consists of an integrated dual one-way valve that guarantees protection of the patient from fluid ingress and simplified user control of the device, allowing: drainage, controlled manual aspiration (pull & push manoeuvre), pleural lavage and talc slurry (micronised talcum powder in sterile injectable solution, saline solution), all without the use of additional components resulting in improved patient management.

## 4. INSTRUCTIONS FOR USE

#### PREPARATION FOR USE

Before using the set, carry out the preparation procedure as described below:

- Check the intact state of the packaging to ensure functionality and sterility of the product.
- Open the package using sterile procedure
- Determine the position for needle insertion with the aid of ultrasound or radiological investigation.
- Prepare a sterile field around the chosen area, taking normal precautions for skin disinfection.
- Remove the components necessary for use from the packaging, such as the main body of the device with the attached catheter; the catheter is also supplied assembled to the body in the version with butterfly connections.
- The Basic version only features the Verres needle. Connect the main body of the device to the rear needle lock connector and close the protective caps (Fig.4).
- In the "pig-tail" version, slide the rigid tube (Fig. 3) on the catheter, towards the end of the catheter until it is completely straightened.
- Insert the needle through the rear access (1) until the handle is fully retracted. The needle and its protective mandrel should protrude a few millimetres from the end of the catheter.
- In the pig-tail version, remove the rigid tube from the catheter (Fig. 3a).
- In the Forty version Insert the needle into the catheter, through the hole indicated by the arrows (Fig. 2), until the handle comes to a stop. The arrow on the handle must come to match that on the catheter and the angle of the handle itself must be directed to follow the curvature of the tube
- Perform local anaesthesia on the site selected for needle introduction using the usual anaesthetics.
- If the doctor considers it necessary, in order to facilitate the introduction of the catheter, make a small incision in the skin.

#### NEEDLE INTRODUCTION AND CATHETER PLACEMENT

After carrying out the preparation described above,

- Hold the device firmly, using both hands; The hand located on the main body is responsible for providing the insertion thrust, directing the needle and keeping the needle handle in contact with the body of the device. The other hand is recommended to be placed close to



Fig. 8a




Fig. 8b

Fig. 9

**Fig. 5a****Fig. 5b****Fig. 6****Fig. 7**

the distal end; grasping the catheter with the thumb and forefinger will conveniently avoid inserting the needle too far into the patient's pleural cavity once it has passed the skin.

- Place the needle tip against the patient's skin and apply gentle pressure until the green visual indicator disappears (Fig.5a).
- Proceed with the introduction very slowly until the green visual indicator reappears; this means that the end of the needle has entered the pleural space and the protective mandrel has returned to its extended position, fulfilling its safety function (Fig. 5b).
- In the Basic version, which is not equipped with a catheter, insertion of the needle completes the positioning procedure.
- In all other versions, move the handle of the device so as to orientate the direction of catheter introduction. To verify that the introduction has taken place at the desired point, i.e. at the area where the fluid to be drained is accumulated, exploratory sampling can be carried out. To do this, connect the syringe to the luer-lock connector on the needle by removing the corresponding cap, retract the piston and visually observe (Fig.6) the entry of liquid into the syringe. If so, proceed with removal of the needle and placement of the catheter, as shown in (Fig.7).
- Begin removal of the needle and, at the same time, complete catheter placement, as shown in Fig.7. In the 'pig-tail' version, removal of the needle will allow the catheter to resume its rounded profile thanks to the shape memory of the material.

#### FLUID EXTRACTION: (PULL&PUSH)

Once the catheter has been positioned, suction of the fluid to be drained can proceed with the aid of a syringe (9) and collection bag (10). To connect these accessories, open the corresponding plugs (2) located on the two luer-lock connections on the main body. Connect the syringe to one of the two connectors and the collection bag to the other.

The choice of which connector to use to connect the syringe or bag is left up to the physician, who will decide on the basis of ease of operation.

Pull the syringe plunger (Fig.8a): fluid will be sucked from the chest cavity into the syringe. By pressing the plunger, the liquid will automatically be sent to the collection bag (Fig.8b). Repeating the pulling and pressing operation in sequence will result in the evacuation of the liquid, which will be collected in the bag.

No additional operation is necessary and the selector switch must remain in the initial position (the knob is oriented according to the axis of the catheter as in Fig. 8 a-b).

At any time, the bag can be replaced or the syringe or both can be removed, depending on the outcome of the drainage. It is advisable, when removing the accessories, to close the protective caps to avoid possible contamination. In any case, if the caps remain open, there is no risk of reflux towards the patient.

In the Basic version, once the needle is positioned, it is possible to proceed with the aspiration of the liquid with the help of a syringe and collection bag. To connect these accessories, open the corresponding plugs located on the two luer-lock connections on the main body. After fluid collection, the needle must be removed as it should never be left in place for a long time and without constant supervision by the doctor.

#### INSERTING THE EXTENSION LINE (for the respective versions)

The device is equipped with a butterfly connection designed to interpose an extension line between the catheter and the main body. In this way, the drain can be easily fixed and kept in place for a long time with greater comfort for the patient. To insert the extension line, the catheter must be clamped using the disposable klemmer (11) included in the package. Then the butterfly is disconnected by connecting the line as shown in Fig. 9. At the end of the operation, the klemmer can be opened and then the catheter can be attached and the fluid collection accessories connected.

#### GRAVITY DRAINAGE: COLLECTION BAG

After the catheter has been positioned and secured, it is possible to proceed with the drainage of fluids by gravity. In this case, the drainage bag included in the package must be connected to any of the available fittings, through which the drained liquids can be collected by gravity.

To drain the fluids, proceed as follows:

- Connect the collection bag to one of the luer-lock fittings (Fig.10a). Before carrying out subsequent operations, make sure that the other connection is closed by the appropriate plug.
- Position the collection system conveniently below the level of the patient.
- Turn the selector switch, located in the main body, as far as it will go, corresponding to a 90° angle to the initial position (Fig.10b). At this point, there will be a direct connection between the catheter and the collection system without the interposition of any valve.
- If you wish to replace the collection bag, e.g. in case of complete filling, turn the selector switch to the home position and then replace the bag. After replacement, rotate the selector again by 90°.

#### FLUID DRAINAGE: POST-SURGERY COLLECTION SYSTEM

After the catheter has been positioned and secured, it is possible to proceed with fluid collection by gravity or suction using an appropriate postoperative collection system (e.g. silicone reservoir, chest drainage system, etc.). In this case, the collection system must be connected to any of the available fittings, possibly using the appropriate connector included in the package (12).

To connect a postoperative drainage device, proceed as follows:

- If a chest drainage system is to be used, connect the connector included in the package to one of the luer-lock fittings (Fig.11).
- Prepare the collection system by carrying out the operations prescribed by the manufacturer, refer to the instructions for use of the respective products; then connect the system to the free end of the connector.
- Turn the selector switch, located in the main body, as far as it will go, corresponding to a 90° angle to the initial position. At this point, there will be a direct connection between the catheter and the collection system without the interposition of any valve.
- If you wish to replace the collection system, e.g. in case of complete filling, turn the selector switch to the home position and then replace the bag or drain. At the end of the operation turn the selector switch by 90° back to the previous position.

#### INFUSION OF SOLUTIONS (PLEURAL TALC SLURRY OR LAVAGE)

After positioning and securing the catheter, it can be used as an access for the infusion of pleural talc solution (micronised talcum powder in a sterile injectable solution) or washing solution (e.g. saline).

To perform this manoeuvre, proceed as follows:

- connect the syringe containing the liquid to be infused to one of the luer-lock connections and ensure that the other connection is closed with the cap (Fig.12).
- Turn the selector switch, located in the main body, as far as it will go, corresponding to a 90° angle to the initial position. At this point, there will be a direct connection between the catheter and the syringe without the interposition of any valve and the infusion of the solution to the patient can be initiated.
- When finished, turn the selector switch to the home position and disconnect the syringe. If there is a need to continue drainage, the operations described above (pull-push, gravity drainage and suction) can be performed; please refer to the previous sections for details.

#### CATHETER FIXING

Inside the set you will find a useful accessory for fixing the catheter (13). You can take advantage of this accessory to secure the placement of the drain, fixing it to the skin by means of a suture.

The fixing system can be positioned close to the percutaneous access, so that it covers a portion of the tube (catheter and/or extension line); one of the two side wings is equipped with a handle and can be tied to the opposite wing by passing underneath the tube: this manoeuvre prevents the tube from dislodging (Fig. 13).

#### DEVICE REMOVAL

At the end of the drainage operation, the device can be removed. For removal, proceed as follows:

- Remove any bandages and dressings and disinfect the area around the point of entry of the drainage tube.
- Remove any fixing suture and prepare the wound closure suture (e.g. tobacco pouch) in advance.
- Pull the drain tube out by exerting moderate traction. When extraction is complete, check the intact state of the tube and then close the wound.

#### 5. KIT COMPOSITION

The kit consists of the following components:

- polyurethane catheter with integrated valve system (Fig.1a-b)
- Verres dilator needle (7)
- 2000 ml collection bag (10), equipped with a valve for air evacuation (14) and a tap for emptying (15);
- one syringe (9);
- one disposable scalpel (16);
- one connection for drainage systems (12);
- one extension line (6);
- disposable klemmer (11);
- catheter fixing set (13);

#### 6. COMPATIBLE ACCESSORIES

The compatible accessories produced by Redax are:

- Catheter fixing system
- Unique fitting for Drentech
- Unique extension line

For codes, please contact [info@redax.it](mailto:info@redax.it).

#### 7. COMPATIBLE DEVICES

The 'UNICO' system is compatible with the following devices manufactured by Redax:

- 2000 mL drain bag
- P-eggy drain thoracic drain valve

For codes, please contact [info@redax.it](mailto:info@redax.it).

## 8. DISPOSAL

When disposing of the device after use, appropriate precautions must be taken and the current legal provisions concerning biologically hazardous waste must be observed.

## 9. MATERIALS USED

Polyurethane (PU) for the drain part in contact with the patient.  
LATEX-AND PHTHALATE-FREE MEDICAL DEVICE  
PVC-FREE PACKAGING

## 10. LABELLING AND SYMBOLS USED

In accordance with ISO standard 15223-1.

Date of last version:  
see last page: (REV: XX-XXXX)

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