REDAX®

IT

VALVOLA DI HEIMLICH

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DESCRIZIONE

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EN

DESCRIPTION

The Heimlich non-return valve is a one-way drainage device which prevents the air reflux towards the patient's chest.

Two versions of the device are available:

Standard Heimlich valve (fig. 1)

This device is composed of the following parts:

1. A flexible body (1) fitted with universal tapered connectors (2) at both ends, so as to be connected to drainage tubes and/or collection devices.
2. A lip valve in flexible material (3), connected to one of the two connectors.
3. An additional connector (4) for the connection to the drainage tubes provided with female luer-lock connector.

Heimlich valve with bag (fig. 2)

This device is composed of the following parts:

1. A standard Heimlich valve (as described in the previous paragraph), the end of which is connected to a collection bag (5), by means rubber connection tube (6).
2. The collection bag has 2000 ml capacity and it is provided with an inlet one-way valve (7).
3. Furthermore, the bag is provided with a filter (8) that allows to evacuate air from the patient outwards.

RECOMMENDATIONS FOR USE

The Heimlich valve is primarily indicated in drainage of the chest cavity by gravity in order to allow the evacuation of liquids and/or air, as well as preventing the backflow towards the patient. This way, the normal gradients of intrathoracic pressure can be restored and it can facilitate a complete pulmonary expansion, by restoring the normal respiratory dynamics. Both versions, with and without bag, are especially indicated in emergency conditions.

Note: This device can only be used in safe conditions in the field of applications and using the procedures indicated on this operating instruction sheet, for the specific type of product. The manufacturer declines all liability relating to improper use or use different to that indicated.

PREPARING THE SYSTEM FOR USE

Perform as follows.

1. Pull out the valve from its sterile wrapping.
2. Check the direction of the arrow marked on the valve: the arrow must be ALWAYS oriented in the direction of the drainage flow, that is, from the drainage tube to the collection system (see fig.1). Regarding the version with bag, the positioning is already correctly preset.
3. Connect the drainage tube to the corresponding end, making sure that the connection is suitably strong.
4. In the standard version, connect the second end to a collection system suitable for the application. The version with bag is already equipped with a 2000 ml capacity collection bag.
5. Place the system, previously prepared as above described, below the level of the wound so as to start the drainage by gravity immediately.

DRAINAGE OF THE PATIENT

Both versions of Heimlich valve allow the drainage by gravity of air and liquids from the patient. The one-way valve (3) prevents the backflow of liquids towards the chest.

- VERSION WITH BAG

The combination of the two one-way valves (3) and (7) and the flexible body (1), allows to create a negativity that can turn out to be useful in order to activate the drainage or to remove potential occlusions from the tube. In order to realize this operation, it is enough to compress and release the valve repeatedly as shown in figure 3. In the event of accumulation of air in the collection bag (5), a filter (8) allows its evacuation as well as prevents, at the same time, the accidental outflow of the liquid. In case of need, it is possible to increase the speed of such evacuation by compressing the bag gently.

CONTRAINDICATIONS

No known contraindications exist for drainage use.

GENERAL PRECAUTIONS

- The product is sterile unless the package has been opened or damaged. Do not use if the package has been opened or damaged.
- Exposure to elevated temperatures and ultraviolet light should be avoided during storage.
- For discarding and disposal of the device, adopt suitable precautions and act in accordance with the provisions of the law in force for biologically hazardous waste.
- Single-use device. Reuse may lead to alteration of performance and risks of cross-contamination. Dispose of the device immediately after use. Do not re-use.

STERILE Ethylene oxide sterilized

MATERIALS USED

PS, PVC, PE, S-E/B-S (version with bag).

LATEX-FREE MEDICAL DEVICE.

WARNINGS/GENERAL PRECAUTIONS

1. Always keep the valve and the bag below the level of the wound so as to achieve the drainage by gravity.
2. Make sure that the valve is connected correctly to the drainage tube, by following the orientation of the arrow marked on the device.
3. Always make sure that the connections are suitably strong .
4. Avoid squashing or bending the connection tube between the valve and the collection bag.

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