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THORACIC CATHETER TROCAR THORACIC CATHETER

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DESCRIPTION

The thoracic catheters are sterile PVC drain tubes fitted with an atraumatic tip and closed terminal connector to facilitate positioning. The Trocar thoracic catheters are sterile PVC drain tubes fitted with an internal stylet that allows positioning using a percutaneous technique. The Trocar catheters may be with a closed or sharp tip to satisfy the various needs and practices of the physician. All the devices are available in various sizes to adapt to the patient's physiology.

INDICATIONS FOR USE

The drains are intended for short-term use in clinical conditions that require evacuation of air and/or liquids that have built up in natural and/or neofornate cavities following surgery. They are indicated for use in association with gravity or suction drainage systems.

SUGGESTIONS FOR USE

Positioning the drain

1. Before positioning the drain, carefully select the type, the position and the required size. These parameters depend on the type of operation, the physiology of the patient and the experience of the surgeon.
2. Using an aseptic technique, carefully take the sterile chest drain tube out of the protective wrapping.
3. Position the drain in the wound following an aseptic procedure based on the protocols and practices in use at the hospital. In the case of the Trocar drain, use the stylet provided for positioning.
4. Use the depth markers indicated on all the drains for correct positioning. In the case of the Trocar drain, remove the stylet after positioning.
5. Secure the drain with a skin stitch and/or a dressing.
6. Connect the positioned tube to the chosen drainage system according to the procedures indicated in the relative instruction leaflet. If necessary, the connector included in the package can be

used.

Removing the drain

1. Remove the tape or stitch or whatever was used to secure the drain.
2. Disconnect the drainage system.
3. Remove the drain by pulling it gently and immediately apply a dressing on the wound.

WARNINGS/PRECAUTIONS

1. **Be extremely careful when positioning the device. The procedure is to be carried out exclusively by an expert physician, especially for the Trocar drains where inexperience may cause serious harm to the patient.**
2. **Do not in any way handle the device or make extra holes with the aid of cutting tools.**
3. **If securing the drain with a stitch, be extremely careful not to suture the tube or cut it with needles or sharp tools. Do not over-tighten the suture as this might cause the device to break.**
4. **After positioning the device, check that the perforated section is fully inserted in the wound and that the connections are perfectly sealed.**
5. **Manually remove the device pulling it gently avoiding abrupt manoeuvres. Do not use metal instruments, such as forceps or devices that may cause the device to break.**
6. **If the device is left in place for very long periods of time, it may be difficult to remove. Be extremely careful during the removal operations.**
7. **Once the device has been removed, check that it is fully intact, an essential condition to exclude that device fragments have been left in the cavity caused by accidental rupture.**
8. **It is recommended not to use the device for more than 29 days.**

CONTRAINDICATIONS

The physician has to assess whether there are any clinical conditions that would make drain positioning inadvisable (e.g. need for an immediate thoracotomy, severe haemorrhage, bacterial contamination of the site, etc.).

GENERAL WARNINGS

The product is sterile if the package is intact. Do not use if the package is open or damaged. Avoid exposure to high temperatures and ultraviolet rays during storage.

These devices can be used in safe conditions only in the field of application and according to the procedures indicated in this instruction leaflet in correspondence to the type of product. The manufacturer declines all responsibility for improper use or different from that indicated.

DISPOSAL

Adopt adequate precautions for elimination and disposal of the device and comply with the provisions of the laws in force on biologically hazardous waste.

SINGLE-USE DEVICE.

Reuse may lead to alterations in performance and cross-contamination risks.

STERILE - Sterilized by ethylene oxide

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