REDAX®

IT

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CONTROLLED-SUCTION DRAINAGE SYSTEM

INSTRUCTIONS FOR USE

DESCRIPTION

The controlled-suction drainage system illustrated in Fig. 1 is available with bellows in different capacities. The suction system is equipped with a dilatory spring (4) located inside the collection bellows (1). The spring limits the dimensions and allows a constant suction level (about 100 mmHg) and collection capacity, guaranteed to the nominal value marked on the bellows. The graduation on the lower part of the system allows checking the quantity of fluid collected. Each system is equipped with an extension tube and relevant multi-gauge connector adaptable to drainage tubes of between 06 and 18 CH diameter. The drainage system can be equipped with a PVC or silicone drainage also available with pre-assembled guide needle.

WARNINGS / PRECAUTIONS

- Avoid loops or squashing the patient drainage tube or the suction tube.
- Do not use if package has been opened or damaged.
- The quantity of fluid collected in the system corresponds to the

indication shown on the graduation only when the system has exhausted its suction action i.e. when it is fully extended.

Note: This device can only be used in safety conditions in the field of applications and using the procedures indicated on this operating instruction sheet for the specific type of product. The manufacturer declines all liability relating to improper use or use different to that indicated.

- Single-use product. Reuse may lead to alteration of performance and risks of cross-contamination.

USING THE SYSTEM

1. Insert the drainage tube (9) into the wound from the outside to the inside using the guide needle of corresponding diameter. Always maintain rigorous asepsis. Fix the drainage tube to the outside of the wound with a suture or plaster.
2. Cut the multi-gauge connector (7) in the position corresponding to the diameter of the drainage tube used (Fig. 1) and insert the distal end of the drainage in the connector.
3. To use two drainages the Y-connector can be used (Fig. 5): Cut the multi-gauge connector (7) at the end of the suction tube (Fig. 5), then proceed for both as in point 2.
4. Screw the connector (3) at the end of the suction tube to the special connector on the cap of the suction system.
5. Close the clamp (5) on the suction tube, check that the trigger valve(2) on the cap of the suction system is open and compress the bellows with both hands (Fig. 3) until full compression is obtained. When the bellows is released the trigger valve (2) closes automatically. If the bellows compression operation has not been carried out in a satisfactory manner, it can be repeated by opening the valve and again compressing the system. When compression has been completed, the valve must be closed completely and will not be used again.
6. Open the clamp (5) on the suction tube to activate suction. A non-return valve prevents the drained fluid from flowing back to the patient. The valve guarantees the passage of fluid from the patient to the bellows preventing the accumulation of clots. On request a version without non-return valve is available.
7. The suction system is equipped with a fixing band (8) for attachment.

SYSTEM REPLACEMENT

When the suction system has exhausted its suction action (when in fully extended position) it has used up its collection capacity and can be replaced. To carry out this operation, proceed as follows:

- Prepare a new suction system.
- Close the clamp (5) on the suction tube of the system in use.
- Unscrew the suction system in use from the relevant suction tube.
- Replace the system with the new one, screwing the suction tube to its Luer-lock connector.
- Follow steps 5 and 6.
- Dispose of the replaced system. For discarding and disposal of the device, adopt suitable precautions and act in accordance with the provisions of the law in force for biologically hazardous waste.

CONTRAINDICATIONS

There are no known contraindications.

STERILE - Ethylene oxide gas sterilized.

The product is sterile unless the package has been opened or damaged. Discard after each single use, do not re-use. Avoid exposure to high temperatures and ultraviolet rays during storage. For discarding and disposal of the device, adopt suitable precautions and act in accordance with the provisions of the law in force for biologically hazardous waste.

Note: The liquid contained in the drainage system must absolutely not be used for infusion.

MATERIALS USED

PEHD, PVC, PP, PS
LATEX-FREE MEDICAL DEVICE. PVC-FREE PACKAGING.

CLOSED-CIRCUIT CONTROLLED-SUCTION DRAINAGE SYSTEM

INSTRUCTIONS FOR USE

DESCRIPTION

The closed-circuit controlled-suction drainage system, illustrated in Fig. 2, is available with bellows in different capacities and is always equipped with a suction tube with relevant multi-gauge connector for a PVC drainage tube (9). The system has the same features as the standard system. The closed-circuit system is fitted with a connector in

connettor in the lower part, to which the mesh is connected for collection (11) with an anti-reflux valve. The mesh can be of different volume. The mesh is used for emptying "garmoški", when it produces its own resource volume, in this way it is possible to restore and prolong suction. It is possible to replace the mesh or empty it by means of the drainage crane (12). The system of the closed type can be equipped with drains from PVC or silicone, and there is also a configuration with a pre-assembled three-way connector.

PRELIMINARY MEASURES PRECAUTIONS

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the lower part to which a collection bag (11) fitted with an anti-reflux valve is connected. The bag may have different capacities. The bag is used to empty out the contents of the bellows when it has used up its collection capacity, so that the suction action can be restored and use prolonged. The bag can be replaced individually or emptied through the drainage tap (12). The closed-circuit set can be equipped with PVC or silicon drainage also available with pre-assembled guide needle.

WARNINGS / PRECAUTIONS

- Avoid loops or squashing the patient drainage tube or the suction tube.
 - Do not use if package has been opened or damaged.
 - The quantity of fluid collected in the bellows corresponds to the indication shown on the graduation only when the system has exhausted its suction action i.e. when it is fully extended.
- Note:** This device can only be used in safety conditions in the field of applications and using the procedures indicated on this operating instruction sheet for the specific type of product. The manufacturer declines all liability relating to improper use or use different to that indicated.
- Single-use product. Reuse may lead to alteration of performance and risks of cross-contamination.

USING THE SYSTEM

To fit the drainage, the connection to the drainage, the Y-connector and the connection to the suction system, see the corresponding points 1, 2, 3 and 4 of the description of the standard controlled-suction drainage system.

5. Follow the instructions for use of the standard bellows taking the precaution of closing the clamp (10) on the bag connection.
6. The suction system is equipped with a fixing band (8) for attachment.

RESTORING THE BELLOWES

When the bellows has used up its collection capacity the fluid can be emptied out into the bag and the suction action of the system restored. This operation can however be carried out at any time, for example, to check the quantity of fluid drained observing the graduation on the bag. To empty the bellows proceed as follows:

- Close the clamp (5) on the suction tube.
- Open the clamp (10) on the bag connection.
- Compress the bellows with both hands (Fig. 3) until it is completely empty.
- When grip is released, the bellows is automatically restored.
- Close the clamp (10) on the bag connection.
- Open the clamp (5) on the suction tube to activate suction.

REPLACING THE BAG

The bag (11) can be replaced at any time. To replace the bag, proceed as follows:

- Ensure that the clamp (10) on the bag connection is closed.
- Disconnect the bag from the bellows by unscrewing the connector.
- Close the bag with its cap.
- Fit a new bag screwing the Luer-lock connector to the bellows.
- The full bag must be disposed of. For discarding and disposal, adopt suitable precautions and act in accordance with the provisions of the law in force for biologically hazardous waste.

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The product is sterile unless the package has been opened or damaged. Discard after each single use, do not re-use. Avoid exposure to high temperatures and ultraviolet rays during storage for discarding and disposal of the device, adopt suitable precautions and act in accordance with the provisions of the law in force for biologically hazardous waste.

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LEGEND

- 01 Collection bellows
- 02 Trigger valve
- 03 Suction tube connector
- 04 Dilatory spring
- 05 Suction clamp
- 06 Suction tube
- 07 Multi-gage connector (06-18 CH)
- 08 Fixing band
- 09 Drainage tube
- 10 Bag clamp
- 11 Collection bag
- 12 Drainage tap

Date of issue of the last version :
see last page : (REV.: XX-XXXX)

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