

P-Eggy drain

Chest Drain Valve

Small and portable
drainage system
providing quick
air leak confirmation



REDAX™

MEDICAL DEVICES

An innovative device that simplifies spontaneous pneumothorax management, and it is the ideal solution for home care treatment. P-Eggy Drain can be directly connected to any chest catheters, providing a safe and cost-effective outcome.

Why do you need P-Eggy Drain?

In case of any air leakage, P-Eggy Drain boosts the drainage effect, and at the same time, prevents any air return to the Patient. Thanks to the integrated well, air leakage may be detected, and the Patient's conditions are under control.

What is P-Eggy Chest Drain Valve?

- 100 ml capacity collection chamber
- Integrated non-return valve provides excellent protection
- Air leakage can be detected through the dedicated air leak well
- Silicone bulb can be squeezed, creating a negative pressure
- Provided with multiple chest tube connectors
- Needleless access sampling port for fluid removal
- A lightweight solution for an easy patient mobility

connectors



air leak well



P-Eggy Drain Chest Drain Valve _____ code **10808**

Packaging: 10 units for box

Connectors kit _____ code **10543**

Packaging: 10 units for box



REDAX™

REDAX s.p.A. Via Galileo Galilei, 18
46025 Poggio Rusco (MN) Italy
Tel. +39 0386 830582 Fax +39 0386 51898
E-mail: info@redax.it - www.redax.it

video tutorial

